# CereForm<sup>®</sup>

**FR** - Implant mammaire pré-rempli de gel de silicone • **EN** - Silicone gel-filled breast implant • **DE** - Mit Silikongel vorgefülltes Brustimplantat • **NL** - Met siliconengel gevuld borstimplantaat • **ES** - Implante mamario relleno de gel de silicona • **PT** - Implante mamário pré-cheio com gel de silicone • **IT** - Impianto mammario con silicone in gel • **EL** - Μαστικό εμφύτευμα προγεμισμένο με τζελ σιλικόνης • **SV** - Bröstimplantat förfyllt med silikongel • **IS** - Gelfylltir sílkon brjóstapúða • **BG** - Гръден имплант с пълнеж от силиконов гел • **HU** - Szilikongéllal töltött emlőimplantátum • **RU** - Имплантат молочной железы с силиконовым гелем • **ET** - Silikoongeeliga eeltäidetud rinnaimplantaat • **LV** - Ar silikona gēlu pildīts krūts implants • **LT** - Krūtų implantas, užpildytas silikono geliu • **TR** - Silikon jel dolu meme implantı • **RO** - Implant mamar umplut cu gel de silicon • **HE** - שתל שד עשוי סיליקון ממולא ג'ל • **AR** - حشوة ثديي مملوءة مسبقاً بهلام السيليكون

|                                            |     |
|--------------------------------------------|-----|
| Français - Notice d'utilisation.....       | 08  |
| English - Instructions for use.....        | 17  |
| Deutsch - Benutzermerkbblatt.....          | 26  |
| Nederlands - Gebruiksaanwijzing.....       | 35  |
| Español - Instrucciones de uso.....        | 44  |
| Portugues - Informações de utilização..... | 53  |
| Italiano - Istruzioni per l'uso.....       | 62  |
| Σλληνικά - Οδηγισο χρρσησ.....             | 71  |
| Svenska - Bruksanvisning.....              | 80  |
| Íslenzka - Notkunarleiddbeiningar.....     | 89  |
| български - брошура за ползвателя.....     | 97  |
| Magyarul - Felhasználói tájékoztató.....   | 106 |
| Русский - Инструкция по применению.....    | 115 |
| Eesti - Kasutajjuhend.....                 | 124 |
| Latviešu - Lietotāja brošūra.....          | 133 |
| Lietuvių - Naudojimosi instrukcija.....    | 142 |
| Türk - Kullanma talimati.....              | 151 |
| Română - Instrucțiuni de utilizare.....    | 159 |
| עברית - הוראות שימוש.....                  | 168 |
| تعليقات الاستخدام - العربية.....           | 176 |

**FR** - Consulter le manuel d'utilisation  
 • **EN** - See user manual • **DE** - Raadpleeg de gebruiksaanwijzing • **NL** - Bitte die Gebrauchsanleitung lesen • **ES** - Consultar el manual de instrucciones • **PT** - Consultar o manual do utilizador • **IT** - Consultare le istruzioni per l'uso • **EL** - Συμβουλευτείτε το εγχειρίδιο χρήσης • **SV** - Se användarhandboken • **IS** - Sjá notendahandbók • **BG** - Моля, прочетете наръчника за употреба • **HU** - Olvassa el a használati útmutatót • **RU** - См. руководство по эксплуатации • **ET** - Tutvu kasutusjuhendiga • **LV** - Skatīt lietošanas rokasgrāmatu • **LT** - Žr. naudojimo vadovą • **TR** - Kullanma kılavuzuna başvurunuz • **RO** - Consultați manualul de utilizare • **HE** - ראה מדריך למשתמש • **AR** - راجع تعليمات الاستخدام

**FR** - Craint l'humidité • **EN** - Protect from moisture • **DE** - Droog houden • **NL** - Vor Nässe schützen • **ES** - Proteger de la humedad • **PT** - Proteger da humidade • **IT** - Teme l'umidità • **EL** - Να προστατεύεται από την υγρασία • **SV** - Fuktkänslig • **IS** - Verjð gegn raka • **BG** - Да се пази от влага • **HU** - Tartsa távol a nedvességtől • **RU** - Беречь от попадания влаги • **ET** - Ei ole niiskuskindel • **LV** - Sargāt no mitruma • **LT** - Saugoti nuo drėgmės • **TR** - Rutubetten olumsuz etkilenir • **RO** - Protejați de umiditate • **HE** - יש להגן מפני לחות • **AR** - أبقها جافة

**FR** - Limite de température  
 • **EN** - Temperature limit • **DE** - Temperaturlimit • **NL** - Temperaturlimiet • **NL** - Temperaturbegrenzung • **ES** - Limite de temperatura • **PT** - Limite de temperatura • **IT** - Limite di temperatura • **EL** - Όριο θερμοκρασίας • **SV** - Temperaturgräns • **IS** - Takmörkun hitastígs • **BG** - Температурна граница • **HU** - Hőszékelny • **RU** - Предельно допустимая температура • **ET** - Piirtemperatuur • **LV** - Temperatūras ierobežojums • **LT** - Temperatūros ribos • **TR** - Sıcaklık limiti • **RO** - Limita de temperatură • **HE** - מגבלת החד האقصى לدرجة الحرارة • **AR** - الحد الأقصى لدرجة الحرارة

**FR** - Ne pas restériliser • **EN** - Do not resterilise • **DE** - Niet opnieuw steriliseren • **NL** - Nicht erneut sterilisieren • **ES** - No reesterilizar • **PT** - Não voltar a esterilizar • **IT** - Non riesterilizzare • **EL** - Να μην επαναποστεριώνεται • **SV** - Får inte omsteriliseras • **IS** - Má ekki smítsæfa aftur • **BG** - Да не се рестерилизира • **HU** - Újrasterilizálása tilos • **RU** - Повторная стерилизация запрещена • **ET** - Mitte korduvsteriliseerida • **LV** - Nesterilizēt atkārtoti • **LT** - Pakartotinai nesterilizuoti • **TR** - Yeniden sterilize etmeyiniz • **RO** - Nu reesterilizați • **HE** - אין לעקר מחדש • **AR** - لا تُعد تعقيمها

**FR** - Ne pas réutiliser / A usage unique  
 • **EN** - Do not reuse / For single use only • **DE** - Niet opnieuw gebruiken / Voor eenmalig gebruik • **NL** - Nicht zur Wiederverwendung / Für den einmaligen Gebrauch • **ES** - No reutilizar / Un solo uso • **PT** - Não reutilizar/de utilização única • **IT** - Non riutilizzare / Monouso • **EL** - Να μην επαναχρησιμοποιείται / Μίας χρήσης • **SV** - För engångsbruk / Återanvänd inte • **IS** - Má ekki endurnota / Einnota • **BG** - Да не се използва повторно / За еднократна употреба • **HU** - Újrahasznosítása tilos/Kizárólag egyszeri használatra • **RU** - Повторное использование запрещено (для одноразового использования) • **ET** - Ei ole korduvkasutatav / Ühekordseks kasutamiseks • **LV** - Nelietot atkārtoti / Tikai vienreizējai lietošanai • **LT** - Pakartotinai nenaudoti / vienkartiniam naudojimui • **TR** - Yeniden kullanmayınız/Tek kullanımlıktır • **RO** - Nu reutilizați/Pentru o singură utilizare • **HE** - לא לשימוש חוזר / לשימוש אחד בלבד • **AR** - لا تُعد استخدامها

**FR** - Orientation verticale • **EN** - This way up • **DE** - Deze zijde boven • **NL** - Vertikale Ausrichtung • **ES** - Orientación vertical • **PT** - Orientação vertical • **IT** - Orientamento verticale • **EL** - Κάθετος προσανατολισμός • **SV** - Vertikal orientering • **IS** - Þetta snýr upp • **BG** - Вертикална ориентация • **HU** - Tartsa álló helyzetben • **RU** - Вверх • **ET** - Püstine paigutus • **LV** - Šajā virzienā uz augšu • **LT** - Laikyti vertikaliai • **TR** - Dikey konumlandırma • **RO** - Cu partea acastă în sus • **HE** - כיוון זה למעלה • **AR** - اتجاه عمودي

**FR** - Limite de gerbage en nombre  
 • **EN** - Stacking number limit • **DE** - Stapellimiet (aantal) • **NL** - Stapelbegrenzung in Zahlen • **ES** - Límite de apilado para almacenamiento • **PT** - Limite de empilhamento em número • **IT** - Limite di accatastamento in numero • **EL** - Όριο στοίβαξης σε αριθμό • **SV** - Högsta antal att stapla • **IS** - Þetta snýr upp • **BG** - Да не се съхраняват поставени един върху друг • **HU** - Korlátozottan rakható halomba • **RU** - Предельно количеству ярусов в штабеле • **ET** - Vínastamismäär • **LV** - Kraušanas ierobežojums pēc skaita • **LT** - Sandėliavimo aukščio riba • **TR** - Sayı olarak istifleme limiti • **RO** - Limita numărului de stivuire • **HE** - מגבלת חד • **AR** - حد عدد التراص

**FR** - Non pyrogène • **EN** - Non-pyrogenic • **DE** - Niet-pyrogeen • **NL** - Nicht pyrogen • **ES** - No pirógeno • **PT** - Não pirógeno • **IT** - Non pirogeno • **EL** - Μη πυρογενές • **SV** - Icke-pyrogen • **IS** - Ekki hitavaldandi • **BG** - Апирогенен • **HU** - Nem gyúlékony • **RU** - Не содержит пирогенных веществ • **ET** - Mittepyrogeenne • **LV** - Nepirogens • **LT** - Nepirogeninis • **TR** - Pirojen olmayan ürün • **RO** - Nepirogen • **HE** - לא פירוגני • **AR** - غير مولدة للحمى

**FR** - Fragile, manipuler avec précaution • **EN** - Fragile, handle with care • **DE** - Breekbaar, voorzichtig behandelen • **NL** - Zerbrechlich, mit Vorsicht handhaben • **ES** - Frágil, manipular con precaución • **PT** - Frágil, manusear com precaução • **IT** - Fragile, maneggiare con cura • **EL** - Ευθραυστο, απαιτείται προσεκτικός χειρισμός • **SV** - Ömtåligt, hantera med försiktighet • **IS** - Brothætt, meðhöndlið með varúð • **BG** - Крехък продукт, да се обработва внимателно • **HU** - Törékeny, óvatosan kezelje • **RU** - Осторожно: хрупкое! • **ET** - Purunemise oht, ettevaatust käsitlemisel • **LV** - Trausls, apejieties ar to uzmanīgi • **LT** - Elgtis atsargiai • **TR** - Kırılgan ürün, dikkatli bir şekilde taşıyınız • **RO** - Fragil, manipulați cu grijă • **HE** - רגיש, יש לטפל בזהירות • **AR** - هشءة، يلزم التعامل معها بعناية

**FR** - Stérile / Méthode de stérilisation utilisant la chaleur sèche / Système à double barrière stérile • **EN** - Sterile / Sterilised using dry heat / Double sterile barrier system • **DE** - Gestérilisiert mit behulp van droge hitte / Steriles Doppelbarriersystem • **NL** - Steril / mit Dampf sterilisiert / Steriel dubbel barrièresysteem • **ES** - Estéril / Método de esterilización con calor seco / Sistema de doble barrera estéril • **PT** - Estéril / método de esterilização através da utilização de calor seco / Sistema de dupla barreira estéril • **IT** - Sterile / Metodo di sterilizzazione tramite calore secco / Sistema sterile a doppia barriera • **EL** - Στείρο / Μέθοδος αποστείρωσης με ξηρή θερμότητα / Αποστειρωμένο σύστημα διπλού φραγμού • **SV** - Steril / Steriliseringssystem med torr värme / Sterilt dubbelbarriärsystem • **IS** - Smitsæft / smitsæft með þurrhita / Sóðthreinsað tvöfalt hindrunarkerfi • **BG** - Стерилен / Стерилизация със суха топлина / Стерилна двойна бариерна система • **HU** - Steril / Száraz hőkezeléssel sterilizált / Steril kettős akadályrendszer • **RU** - Стерильно: стерилизация сухим теплом / Стерильная система с двойным барьером • **ET** - Sterilne / Steriliseerimismeetod kuiva kuumuse kasutamisel / Steriilne topeltbarjääri süsteem • **LV** - Sterils / Sterilizēts, izmantojot karstu gaisu / Sterila dubultbarjeru sistēma • **LT** - Sterilu / sterilizuota sausu karščiu / Sterili dviguba barjerų sistema • **TR** - Steril/Kuru ısı ile sterilizasyon yöntemi / Steril çift bariyer sistemi • **RO** - Steril/Sterilizat cu ajutorul căldurii uscate/Sistem de barieră dublă sterilă • **HE** - עבר עיקור באמצעות חום • **AR** - معقمة / يتم تعقيمها باستخدام الحرارة • **IS** - מערכת מחום סטרילי ככול الحرارة / نظام حاجز معقم مزدوج

**FR** - Utiliser jusqu'au • **EN** - Use until • **DE** - Houdbaar tot • **NL** - Zu verwenden bis • **ES** - Utilizar hasta el • **PT** - Utilizar até • **IT** - Utilizzare fino al • **EL** - Χρήση έως • **SV** - Använd fram till • **IS** - Notið til • **BG** - Краен срок на употреба • **HU** - Minőségét megőrzi • **RU** - Годно до • **ET** - Kõlblik kuni • **LV** - Izmantot līdz • **LT** - Naudojimo terminas • **TR** - Son kullanma tarihi • **RO** - A se utiliza până la • **HE** - תאריך תחולת שימוש • **AR** - للإستعمال لغاية

**FR** - Attention / Précautions • **EN** - Caution / Precautions • **DE** - Vorichtig • **NL** - Achtung / Vorsichtsmaßnahmen • **ES** - Atención / Precauciones • **PT** - Atenção / Precauções • **IT** - Attenzione / Precauzioni • **EL** - Προσοχή / Προφυλάξεις • **SV** - Varning / Försiktighetsåtgärder • **IS** - Varúð / Varúðar ráðstafanir • **BG** - Внимание / Предпазни мерки • **HU** - Figyelem / Óvintézkedések • **RU** - Внимание! / Предостережения • **ET** - Tähelepanu / Ettevaatusabinõud • **LV** - Uzmanību! / Piesardzības pasākumi • **LT** - Dėmesio / Atsargiai • **TR** - Dikkat / Önemler • **RO** - Atenție/Precauții • **HE** - אמצעי • **AR** - تحذيرات / احتياطات • **IS** - Zehirat

**FR** - Conserver à l'abri de la lumière du soleil • **EN** - Store away from sunlight • **DE** - Niet blootstellen aan zonlicht • **NL** - Vor Sonnenlicht schützen • **ES** - Conservar al abrigo de la luz del sol • **PT** - Conservar afastado da luz do sol • **IT** - Conservare al riparo dalla luce del sole • **EL** - Φυλάσσετε το προϊόν μακριά από την ηλιακή ακτινοβολία • **SV** - Förvara borta från solljus • **IS** - Geymið fjarri sólarljósi • **BG** - Да се съхранява на място, защитено от слънчева светлина • **HU** - Napsugárástól védve tárolandó • **RU** - Хранить в темном месте. Беречь от попадания солнечных лучей. • **ET** - Hoiustamine kaitstuna päikesevalguse eest • **LV** - Neuzglabāt tiešos saules staros • **LT** - Saugoti nuo saulės spindulių • **TR** - Güneş ışınlarından koruyunuz • **RO** - A se feri de lumina soarelui • **HE** - לאחסן מהחום • **AR** - تخزين بعيداً عن أشعة الشمس

**FR** - Ne pas utiliser si l'emballage est endommagé • **EN** - Do not use if packaging is damaged • **DE** - Niet gebruiken wanneer de verpakking beschadigd is • **NL** - Bei beschädigter Verpackung nicht verwenden • **ES** - No utilizar si el envase está dañado • **PT** - Não utilizar se a embalagem estiver danificada • **IT** - Non utilizzare se l'imballaggio è danneggiato • **EL** - Μην χρησιμοποιείτε το προϊόν εάν η συσκευασία έχει υποστεί ζημιά • **SV** - Använd inte om förpackningen är skadad • **IS** - Notið ekki ef umbúðir eru skemmdar • **BG** - Не използвайте, ако опаковката е повредена • **HU** - Ha a csomagolása sérült, ne használja • **RU** - Не использовать, если упаковка продукта повреждена • **ET** - Mitte kasutada, kui pakend on kahjustatud • **LV** - Neizmantoj, ja iepakojums ir bojāts • **LT** - Nenaudoti, jei pažeista pakuotė • **TR** - Ambalajı hasar görmüş ise kullanmayınız • **RO** - Nu utilizați dacă ambalajul este deteriorat • **HE** - אל להשתמש אם לא נשדמה • **AR** - لا تستخدمها في حالة العبوة

SN

FR - Numéro de série • EN - Serial number • DE - Seriennummer • NL - Seriennummer • ES - Número de serie • PT - Número de série • IT - Numero di serie • EL - Αριθμός σειράς • SV - Serienummer • IS - Raðnúmer • BG - Сериен номер • HU - Sorozatszám • RU - Серийный номер • ET - Seerianumber • LV - Sērijas numurs • LT - Serijos nr. • TR - Seri numarasi • RO - Număr de fabricație • HE - מספר סידורי • AR - الرقم التسلسلي -

LOT

FR - Numéro de lot • EN - Batch number • DE - Partijcode • NL - Chargenbezeichnung • ES - Número de lote • PT - Número de lote • IT - Numero di lotto • EL - Αριθμός παρτίδας • SV - Partinummer • IS - Lotunúmer • BG - Партиден номер • HU - Tételszám • RU - Номер партии • ET - Partii numbe • LV - Partijas numurs • LT - Partijos nr. • TR - Parti numarasi • RO - Număr de lot • HE - מספר אצווה • AR - رقم الحزمة

UDI

FR - Code UDI • EN - UDI code • DE - UDI-Code • NL - UDI-Code • ES - Código UDI • PT - Código UDI • IT - Codice UDI • EL - Κωδικός UDI • SV - UDI-kod • IS - UDI-kóði • BG - UDI код • HU - UDI-kód • RU - Код UDI (уникальный идентификатор изделия) • ET - UDI kood • LV - UDI kods • LT - UDI kodas • TR - UDI kodu • RO - Cod UDI • HE - UDI קוד • AR - رمز UDI

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## 1. PRESENTATION OF THE DEVICE

### General information

This leaflet applies to CEREFORM® breast implants manufactured by EUROMI Biosciences and is intended for surgeons performing breast implant procedures. This medical device must be used only by competent, trained and qualified medical practitioners in the normal course of their work.

CEREFORM® breast implants are long-term implantable medical devices pre-filled with silicone gel. These breast implants are sold sterile following dry heat sterilisation. They are for single use and for one patient only. Under no circumstances should they be re-used or re-sterilised.

CEREFORM® breast implants are packed in disposable double packaging to provide a dual microbiological layer and guarantee sterility until the procedure. The cardboard packaging provides additional mechanical protection to the double packaging.

Each CEREFORM® breast implant is traceable through a unique identification number engraved on the patch. CEREFORM® breast implants obtained CE marking in 2021.

### Composition

The raw materials used to manufacture CEREFORM® breast implants belong to the « Polydimethylsiloxane » group, hence CEREFORM® breast implants contain a small amount of « cyclotetrasiloxane - D4 Siloxane », « cyclopentasiloxane - D5 Siloxane » and « cyclohexasiloxane - D6 Siloxane ». All raw materials are medical grade and biocompatible.

CEREFORM® breast implants are made of:

- a **shell** composed of successive layers of silicone elastomers including a barrier layer limiting gel perspiration,
- an **occlusion patch**,
- a **filling gel**.

**NB:** Small amounts of D4, D5, D6 Siloxane and platinum can pass through the shell of a breast implant that is intact. The clinical effects on humans are not yet known but EUROMI Biosciences takes all the necessary steps to measure the potential toxicity of raw materials used to manufacture CEREFORM® breast implants, especially by performing controls to detect any presence of heavy metals, D4, D5, D6 Siloxanes and residual solvents. Platinum is a metal present in the raw materials used to manufacture our breast implants. Therefore, small residual quantities of platinum may penetrate the body, either by diffusing through the intact shell or following an implant rupture. Thanks to the literary review and other available data we can conclude that the platinum contained in the breast implants has a zero oxidation state and presents the lowest toxicity. This platinum does not, therefore, represent a significant risk for patients having breast implants. A brief summary of the main scientific studies carried out with silicone gel-filled breast implants is available on the FDA website ([www.fda.gov](http://www.fda.gov)).

### Description

CEREFORM® breast implants come in a **round shape** that provides volume to the upper part of the breast for a more plunging neckline. It is particularly suitable for women whose breasts are already formed. The round shape is available in several profiles (low, medium, intermediate, high and very high) and the volumes range from 100 to 900 cc.

CEREFORM® breast implants are available with two different surfaces:

- **The smooth surface**, to facilitate implantation and explantation.
- **The micro-texturing**, to improve cell colonisation and reduce the risks of capsular contracture while still facilitating implantation and explantation.

In addition, CEREFORM® breast implants are available with three filling gels: **CEREFORM® Classic**, **CEREFORM® Aptima** and **CEREFORM® Ellipse**.

The CEREFORM® Classic gel is the softest in the range, for a more natural touch, the CEREFORM® Aptima gel is firmer for better implant holding and the CEREFORM® Ellipse gel offers a softer touch and easier implantation thanks to its shape memory.

To choose the most appropriate device for each patient, you should refer to the sales catalogue provided by EUROMI Biosciences.

## Advantages

CEREFORM® implants have several advantages:

- a **highly resistant patch**, undetectable to the touch with a silicone gel barrier layer for maximum safety,
- **dry heat sterilisation**, guaranteed by highly robust ergonomic double packaging designed for easy use,
- **two types of surface** to choose between easier implantation and explantation (smooth surface) and a lower risk of capsular contracture (micro-texturing surface, with no added texturing agent).

## 2. MEDICAL PURPOSE / INDICATIONS

### Medical purpose

In plastic surgery, the aim of a breast implant procedure is to improve the aesthetic appearance of the breast or to increase its volume, while in reconstructive surgery, the aim of the breast implant procedure is to « remodel » the breast.

### Indications

CEREFORM® breast implants are indicated for:

- **breast reconstruction** to restore the breast following a mastectomy (after cancer or trauma). This is a chance for the woman to recover her body image and forget her disease;
- **aesthetic breast augmentation** to increase the breast size and restore the woman's self-esteem, self-confidence and enhances her sense of womanhood;
- **plastic breast augmentation** to correct various congenital or acquired abnormalities such as asymmetry, amastia, aplasia, hypomastia and hypoplasia;
- **replacement** of an old or defective implant.

It is recommended to wait until the end of puberty before having a breast implant procedure, and to comply with the minimum legal age in the country of the procedure, if applicable.

## 3. CONTRAINDICATIONS

CEREFORM® breast implants are contraindicated in the following situations:

- Pre-existing pathology in the implant area
- If the surgeon deems that the patient is in poor physiological condition
- A general infection or infection in the implant area
- Grossly positive axillary involvement / or chest wall involvement
- Progressive breast cancer large tumours (>5 cm), late cancer stage and deep tumours
- Current pregnancy or planned in the near future
- Diabetes
- Known hypersensitivity and / or allergy to silicone
- History or presence of an autoimmune disease
- Psychological instability
- Repeated failure with the implantation of similar implants
- Glycated haemoglobin (HbA1c > 7.5%)
- High risk of cancer recurrence
- Insufficient tissue or fat
- Currently breastfeeding or planned in the near future
- Obesity or a Body Mass Index >40 kg/m<sup>2</sup>
- Smoking
- History or current treatment using radiation with lower-pole scarring and thin, poorly vascularized skin/tissue / microwave diathermy / or steroids.

CEREFORM® breast implants are not suitable for men (including sex-change cases).

## 4. POSSIBLE COMPLICATIONS

Like any surgical procedure, breast implant procedures involve surgical and post-operative risks. The practitioner is responsible for the patient's pre-operative assessment and the surgical method used. He/she must inform the patient of the risks associated with the procedure and potential post-operative complications. The practitioner must obtain the informed consent of the patient by asking her to sign the form before the procedure. He/she must also inform the patient about alternative methods to implant procedures.

### Breastfeeding problems

Following breast augmentation, some women may have difficulty breastfeeding their children. Similarly, women who have undergone a mastectomy and breast reconstruction surgery may not be able to breastfeed due to loss of breast tissue and the glands that produce milk.

### Asymmetry

Post-operative asymmetry is the result of choosing the wrong implant (size, shape) that is disproportionate to the other breast, or a different tissue reaction between the two breasts. If there is considerable asymmetry and the patient is dissatisfied, the implant can be removed or replaced. Asymmetry appearing several months or years after the procedure suggests capsular contracture or a ruptured implant. In this case, a thorough examination is necessary, and the implant may have to be removed.

### Atrophy of the breast tissue / deformation of the chest wall

The pressure exerted by the implant may cause the breast tissue to become thinner and retract, increasing the visibility and palpability of the implant. Atrophy of the breast tissue may potentially deform the chest wall.

### Calcification

Calcification is when calcium is deposited in the tissue around the breast implant. These deposits, although quite rare, are painful and can damage the implant, which must then be explanted.

### Breast cancer

The number of breast cancer cases in women with (a) breast implant(s) is not related to the number of breast cancer cases in the general female population.

### Change in the sensitivity of the breast or nipple

A breast implant procedure can increase or decrease breast and/or nipple sensitivity. This change in sensitivity can be temporary or permanent.

### Delayed healing / hypertrophic scars

When the incision does not heal well this can result in unsightly, hypertrophic or keloid scars. Delayed healing may also occur in the following cases: practice of some sports, an infection, sutures that are too tight, an implant that is too large. If the problem persists, the scars may have to be treated surgically.

### Capsular contracture

The fibrous capsule, which forms naturally around any foreign material implanted in the human body, can retract around the latter, causing abnormal compression. This is painful and can deform the breast and rupture the implant. The implant may have to be removed (with or without reimplantation). It is strongly advised not to perform a capsulotomy to treat capsular contractures as this procedure can cause wrinkles or even rupture the shell.

### Deflation of the implant

Implants deflate as a result of a rupture. This is rare with silicone gel implants because of the gel's cohesion. Any deflation detected by the patient should be interpreted as a ruptured implant and more thorough examinations will be necessary. If it is established that the implant has ruptured, it must be explanted.

### Pain

Post-operative pain due to the surgical procedure can be experienced two to three days after the surgery. The post-operative pain intensity varies according to the patient. The patient then experiences discomfort over the following month. These pains can be treated with analgesics. If there is any persistent or new pain in the area of the implant, the patient should be examined to rule out a complication.

### Effects on children

Although there is currently no established method to accurately detect silicone levels in breast milk, a study measuring silicon levels (one component of silicone) did not indicate higher levels in the breast milk of women having silicone gel-filled implants than in women without an implant.

### Permanent implant explantation

If multiple implant-related complications occur repeatedly, or if the surgeon deems that the patient's condition requires the explantation of the implant, permanent implant explantation should be considered, meaning that the implant cannot be replaced and leading to unsightly results (sagging breasts, wrinkles).

### **Extrusion**

When there is an excessive pressure on the implant, it can protrude through the incision or the skin.

### **Gel leakage / silicone migration**

Despite its barrier effect, the silicone shell containing the silicone gel is not perfectly impermeable. Small amounts of silicone can therefore diffuse out of the implant and spread into the tissues. Silicone gel is not toxic to the body but local reactions can occur, forming small fibrous capsules.

### **Galactorrhoea**

Galactorrhoea is the flow of milk through the nipple when the woman is not breastfeeding. This can occur on both sides (bilateral) or in one nipple only (unilateral).

### **Haematoma**

To prevent a haematoma in the implant area, meticulous haemostasis must be performed during the procedure. In case of persistent haematoma, a puncture can be performed taking all necessary precautions to avoid damaging the implant. Appropriate medical compression in the implant area for a few weeks following the procedure will reduce postoperative oedema.

### **Implant noticeable to the touch / visible implant**

Improper initial implant positioning, movement of the implant, inappropriate an implant volume, an inappropriate implant type and cohesion are all factors that can make the implant noticeable to the touch and/or visible. If the patient is dissatisfied, a new procedure may be considered.

### **Infection / inflammation / fibrosis**

Short-term or long-term post-operative infections are only very rarely reported with breast implant procedures. However, any infection should be treated as soon as it occurs. If antibiotics are unable to treat the infection, the implant can be explanted. Inflammation may manifest as redness, swelling, a sensation of heat or a throbbing pain. Fibrosis occurs following substantial tissue degradation or when inflammation develops in an area where tissues are not regenerated.

### **Lymphoedema or lymphadenopathy**

A breast implant procedure can cause swelling or an abnormal reaction in one or more axillary lymph nodes.

### **Connective tissue disease**

No relation has been established between silicone gel-filled breast implants and connective tissue diseases, breast cancer or reproductive problems but these complications cannot be ruled out as studies on these subjects are not sufficiently relevant.

### **Poor positioning / movement**

The implant may move because of poor initial positioning, trauma in the implant area or significant early slackening of the surrounding tissue, failing to hold the implant sufficiently. This results in loss of functionality of the implant (implant herniation, reversal or a change in breast shape) requiring further surgery.

### **Mass, lump, cyst, granuloma**

A lump or mass of inflammatory cells surrounding a foreign substance may occur due to long-term inflammation.

### **Necrosis**

Tissue necrosis can be caused by:

- soaking the implant in an iodine solution before the procedure,
- an abnormal local tissue reaction possibly due, for example, to infection or treating the tissue with radiotherapy prior to the implant procedure, etc.,
- very taut tissue due to insufficient tissue or an overly large implant.

### **Additional surgery**

In view of the various possible complications related to a breast implant procedure, further surgery is always possible. In addition, the limited lifespan of the implant may result in another procedure to maintain the desired results. The patient must understand and accept the risks of further surgeries before making the decision to have an implant procedure.

### **Breast ptosis**

This phenomenon of sagging breasts normally results from ageing, pregnancy or weight loss. Like natural breasts, a breast fitted with a breast implant may sag over years due to tissue distension in the implant area. Ptosis is not dangerous and can be treated surgically.

### **Delays in breast cancer diagnosis**

Having (a) breast implant(s) can delay the diagnosis of breast cancer. This is because the breast implant hides part of the breast and can therefore hinder the interpretation of mammogram results.

### **Unsatisfactory results (size, shape, appearance)**

After a breast implant procedure, the patient may not be satisfied with the general appearance, depending on the style, shape or size of the implant used.

### **Wrinkles / creases / ripples**

It is possible for the implant shell to wrinkle or undulate, forming ripples depending on how well it is held in the pocket and its position in relation to the pectoral muscle according to the surgical indication. The wrinkles can be noticeable on the surface of the skin. The only solution to correct this is explantation.

### **Risks of the surgery / iatrogenic injury, damage**

A breast implant procedure involves risks inherent to the procedure itself, such as the risks and complications associated with a general anaesthesia, iatrogenic injury or damage. It is vital to take all these risks into account in the pre-operative assessment and to inform the patient before surgery.

### **Redness / ecchymosis / bruising**

Bleeding during surgery can cause a change in skin colour. This is a symptom expected with surgery and is temporary.

### **Deep vein thrombosis.**

Refers to the formation of one or more blood clots.

### **Rupture**

The implant can rupture due to surgical trauma (implant damaged during insertion or by surgical instruments) or post-operative trauma (violent impact, excessive compression in the area of the breasts) or natural ageing. This rupture can be asymptomatic (silent rupture), cause the implant to deflate, or change the shape or the appearance of the breast. If there is any doubt, a diagnostic examination must be carried out (mammogram, ultrasound or MRI) to ensure that the implant is intact. Lastly, regular monitoring enables any rupture to be detected early on. If the rupture is confirmed, the implant must be explanted. To guarantee constant mechanical properties and limit the risk of rupture, regular tests are performed according to current standards.

### **Seroma**

This is the accumulation of lymph around the implant, resulting in a temporary increase in the volume of the breast, which usually resolves itself if it is small, or which can be treated by aspiration if there is a significant increase.

### **Breast Implant Associated ALCL (BIA-ALCL)**

Breast Implant Associated ALCL (BIA-ALCL) is a rare type of lymphocytic lymphoma. According to the global regulatory agencies and medical literature, an association has been identified between breast implants and the development of anaplastic large cell lymphoma (ALCL). Indeed, patients with breast implants have a very small but increased risk of developing Breast-Implant Associated ALCL within the scar capsule adjacent to the implant.

Most confirmed cases of BIA-ALCL have occurred in patients with textured surface implants, although there are known cases in patients with smooth-surface breast implants.

You should consider the possibility of BIA-ALCL when a patient presents with late onset, peri-implant seroma. In some cases, patients presented with capsular contracture or masses adjacent to the breast implant. When testing for BIA-ALCL collect fresh seroma fluid and representative portions of the capsule and send for pathology tests to rule out BIA-ALCL including cytological evaluation of seroma fluid or mass and cell block immunohistochemistry testing for Cluster of Differentiation (CD30) and Anaplastic Lymphoma Kinase (ALK) markers.

If your patient is diagnosed with BIA-ALCL, develop an individualized treatment plan in coordination with the patient's multi-disciplinary care team. Most of the published information about treatment describes removal of the implant and the capsule surrounding the implant, and in some patients, treatment with chemotherapy and radiation.

**Faced with functional or physical signs (effusion, volume increase, pain, inflammation, mass, ulceration, deterioration of the general state) occurring in particular at a distance from the post-operative phase in a woman carrying breast implant, the diagnosis of anaplastic large cell lymphoma associated with a breast implant must be evoked.**

## **5. USE OF THE MEDICAL DEVICE**

### **General information**

CEREFORM® breast implants can only be used by competent practitioners with experience in breast surgery. CEREFORM® breast implants must only be used in the operating theatre in a sterile environment.

## Implant insertion


- Check that the double blister is intact. If there is any doubt as to whether either piece of packaging is in perfect condition, do not use the implant.
- Check that the implant has been stored under the conditions specified on the packaging.
- Check that the expiry date has not been exceeded.
- Handle the implant aseptically / Handle the implant with gloves.
- Handle the implant with care. Any implant that has fallen on the floor or been subjected to an impact should not be used.
- The implant must not come into contact with any other substances (glove talc, fabric, iodine solution, etc.).
- It is possible to immerse the implant in a sterile saline solution bath at body temperature to facilitate insertion.
- The device must not be damaged in any way. Users must never add a filling fluid into the shell.
- Take care to avoid sharp and pointed objects when handling the implant.
- Ensure that there are no wrinkles in the shell once it is positioned.
- The incision size must be adapted to the implant size to avoid overly deforming the implant which could damage it and to avoid damaging tissue due to a lack of space.
- Do not use medicinal substances with the implant.
- The practitioner is advised to have a spare implant available to replace any faulty implant or in case of a handling error during surgery.
- The surgeon is advised to have several breast implant sizes available in the operating theatre to have a certain amount of flexibility in the choice of implants.
- To avoid damaging the implant, do not massage or apply needle treatment (acupuncture) to the implant area after the procedure.
- Should the shell rupture, clean and rinse the pocket thoroughly with sterile saline solution to remove any potential traces of gel.
- Any damaged implant must be explanted.

## Surgical techniques

Several surgical techniques can be used to insert a silicone gel-filled breast implant. The surgeon should evaluate the size, projection and surface of the implant to determine the most appropriate incision area and implant position for the patient's anatomy and expected aesthetic results.

For the insertion of CEREFORM® breast implants, the surgeon can make the incision as follows :

- **The periareolar incision (1)** is made around the areola. The junction between the darker areola and the lighter breast skin hides the adjacent scar. This incision has the advantage of healing particularly well. To consider such an incision, the areola must be wide enough to enable the surgeon to insert the implant comfortably. This technique involves a slightly higher risk of deteriorating the sensitivity of the nipple and areola.
- **The inframammary incision (2)** is the technique most widely used for breast augmentation. Made in the skin fold, it leaves a scar that is concealed under the breast. Among the various possible incisions, this one gives the surgeon the best control during the procedure.
- **The axillary incision (3)** is made in the armpit and has the advantage of being located in an area that is less visible. However, in case of further surgery (to change the implant or following complications), an incision may be required elsewhere. This is also the position that gives the least control during the procedure. The surgeon may have to use specific instruments (endoscope, cold light) for better control.


## Positioning of the implant

CEREFORM® breast implants can be positioned as follows:

- **Retroglandular or prepectoral position (1)** : the implant is positioned immediately behind the mammary gland. This type of insertion is technically simple and enables good results when the gland is thick enough to hide the implants. More natural and physiological, this positioning deepens the inframammary fold as the implant's weight is added to one of the breasts. This breast implant position gives the breasts greater mobility on the bust, creating a more natural appearance. However, the results can change if this thickness decreases over the patient's life (pregnancy, menopause, weight loss). In addition, this position makes mammogram interpretations more complex.

• **Retromuscular or retropectoral position (2)** : this position is suitable for thin patients with very small breasts that do not sag. This position places the muscle between the implant and the exterior, making it less noticeable to the eye and to the touch. In this position, mammograms are easier to interpret. However, when the pectoral muscles contract, this may cause the implants to move outwards.

• **Dual plane position (3)** : the implant is placed behind the muscle in the upper part of the breast and behind the gland in the lower part. This dual position combines the advantages of the two previous techniques. It has the advantage of interposing the muscle naturally on the upper section and of integrating the implant onto the breast in the lower section.

### Explanation

- It is recommended to explant the implant through the same incision initially made for insertion.
- If the implant has ruptured, cleaning the breast pocket with a sterile saline solution is highly recommended. If implant gel is present in the pocket, it must be physically explanted by the surgeon.

### Disposal of the product

Any implant that is removed or non-sterile must be disposed of with biological waste involving an infectious risk.

## 6. WARNINGS TO PATIENTS

Breast implants have a limited lifespan. Given the potential complications mentioned above and the natural wear and tear of the implant within the body (daily mechanical stress), the implant may need to be explanted or replaced, involving further surgery. Given the activities and the many factors that can affect the lifespan of the implants, the lifespan of a breast implant cannot be estimated accurately for any given patient. **The Lifespan of CEREFORM® breast implants is estimated at 75% after 9 years of implantation** (Kaplan Meier method).


**The practitioner must warn his/her patients about any potential complications during and after the procedure and inform them of alternatives to a breast implant procedure** (external breast prostheses, breast reconstruction using their own tissue, etc.).

The patient must be informed that if the implant is explanted and not replaced, the result will be unsightly (sagging breasts, wrinkles, etc.). To perform the procedure, **the practitioner must obtain the free and informed consent of the patient on the appropriate form.**

**The patient must be informed of the following:**

- She must carry the completed patient card given to her by the practitioner at all times to facilitate emergency medical care (for example in case of a road accident).
- She must go to check-up appointments and report any trauma, deflation or pain in the implant area.
- The patient must consult a doctor for routine monitoring in order to detect breast cancer.
- The patient must consult a doctor or surgeon if she suspects a complication, especially if there has been any trauma or compression (for example, caused by a violent breast massage, sport or the use of a seat belt).
- The patient should mention the presence of a breast implant in case of imaging examinations (mammogram, ultrasound, MRI) or a medical or surgical procedure near the implant area.

- The patient should consult a doctor or pharmacist before applying topical medications (such as steroids) to the breasts.
- The patient must be cautious concerning some sports that may cause impacts to the implant (combat sports, etc.) and should be aware of the risks of early rupture in case of an impact.
- The patient should be cautious as to the use of cold treatments (e.g. cryotherapy), which can weaken breast implants.
- The patient must not have acupuncture in the breasts.
- The patient must wait at least 3 months after the procedure before becoming pregnant.
- She must examine her breasts once a month, preferably one week after the end of her period or on a set date if she is post-menopausal. The purpose of this self-examination is to detect a lump or any change in the appearance of her breasts. The self-examination is carried out as follows:
  - With the arm raised, use the tips of the three middle fingers of the free hand to examine the opposite breast.
  - Make small circular movements without removing the fingers from the skin.
  - Examine the whole breast: from the clavicle to the base of the breast and from the sternum to the armpit.
  - Move upwards and downwards repeatedly.

## 7. IMPLANT ASSESSMENT METHOD AND FREQUENCY

The recommended consultation frequency is 1 month, 3 months, 6 months and 12 months after the procedure then every year, even if there are no implant-related symptoms. In case of implant deflation, trauma since the last examination, implant deformation, pain experienced by the patient, or any other signs that may indicate a capsular contracture or a rupture of the implant, a mammogram, ultrasound or MRI should be performed. An imaging examination to detect a rupture is also recommended every year from the eighth year after the procedure.

As silicone is partially radiopaque, the implant can obscure the area underlying the implant during medical imaging (mammogram, ultrasound, MRI). The radiology technician must adapt the technique to choose the most appropriate angles of vision as some areas may be obstructed by the implant. For mammograms, the technique must always be adapted because there is also a risk of the implant rupturing or being weakened.

## 8. TRANSPORTATION AND STORAGE OF THE MEDICAL DEVICE

CEREFORM® breast implants must be transported and stored in their cardboard packaging, away from light, at a temperature between 0°C and 40°C, away from humidity and in a horizontal position. Handle with care.

## 9. RESTERILISATION AND RE-USE

CEREFORM® breast implants are sterile and for single use. CEREFORM® implants must never be re-used. Re-using a non-sterile CEREFORM® implant may result in serious clinical complications and potentially the patient's death. Similarly, CEREFORM® implants must never be re-sterilised as this would greatly deteriorate its mechanical properties and dimensions.

## 10. LABELS, PATIENT CARD AND SURGEON'S STATEMENT

Each implant is supplied with:

- **Eight labels:** these show the essential information about the implant: serial number, volume, texturing and implant type, commercial reference, batch number, brand, location (left/right). The surgeon must add the patient's name and the procedure date to these eight labels. One label must be affixed to the patient card,

one to the statement, and the others must be kept in the patient's medical record.

- **A patient card:** the patient must keep this with her at all times as it identifies the implants she is carrying.
- **An informed consent document:** it must be signed by the patient to prove that she is aware of the risks associated with the procedure and with carrying breast implants.
- **The surgeon's certificate:** this statement declares that the surgeon has informed the patient of complications inherent to the breast implant procedure and that the patient has signed the informed consent form. This statement must be completed by the surgeon and returned to EUROMI Biosciences to guarantee the traceability of the implant.

## 11. MEDICAL DEVICE VIGILANCE CASES

Any serious incident or risk of a serious incident that has resulted or is likely to result in the death or severe deterioration in the health of a patient, user or third party involving a CEREFORM® breast implant must be reported to the appropriate authorities immediately and to EUROMI Biosciences: [materiovigilance@euromi-biosciences.com](mailto:materiovigilance@euromi-biosciences.com).

## 12. HANDLING RETURNS

Explanted products for which a complaint has been filled or that have caused a serious incident or risk of a serious incident should be reported and returned to the EUROMI Biosciences local medical device vigilance representative. Before returning the removed product to the manufacturer, it must first be decontaminated and disinfected (in accordance with the health facility's current procedures). The implant should not be returned to EUROMI Biosciences if the patient is infected with HIV or hepatitis or if she is a known or suspected carrier of another infectious agent.

## 13. GUARANTEES – WARRANTY TERMS

EUROMI Biosciences guarantees that every precaution has been taken in manufacturing CEREFORM® breast implants. Nevertheless, EUROMI Biosciences cannot control the conditions of storage and use of its products, nor the choice of patients or surgical technique. Similarly, EUROMI Biosciences cannot guarantee the outcome of the procedure or the patient's satisfaction and cannot foresee the occurrence of any of the adverse reactions described in this leaflet. As such, EUROMI Biosciences cannot be held liable in case of an incident due to the use of the product or in case of loss of or damage to any of its products. EUROMI Biosciences undertakes to replace an implant if and only if it was defective at the time it left its premises. This clause acts as a warranty and excludes any other warranties that are not stated in the text above, either explicit or implicit under the terms determined by law, or otherwise, including but not limited to any implicit guarantee of marketable quality or suitability for use.

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